


Anti-Integrin Alpha 4 (Itga4/Cd49D) Rabbit Polyclonal Antibody For Wb, Elisa - P13612

Item Number: CM0017999



Introduction

Rabbit polyclonal antibody targeting human Integrin alpha 4 (ITGA4/CD49d), validated for Western blot (WB) and ELISA. Suitable for cell adhesion and immune cell trafficking research.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 4 (ITGA4/CD49d) Rabbit pAb
Remarks/Alias	IA4; CD49D; Integrin alpha 4 (ITGA4/CD49d)
Species	Human
Gene ID (Human)	3676
Gene ID	3676
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 770-970 of human Integrin alpha 4 (ITGA4/CD49d) (NP_000876.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13612
Protein Weight	115kDa
Shipping	Ice bag