

Anti-Integrin Alpha 5 (Itga5/Cd49E) Monoclonal Antibody () For Fc - P08648

Item Number: CM0010106



Introduction

Rabbit monoclonal antibody conjugated to , targeting Integrin alpha-5 (ITGA5/CD49e). Reacts with human and mouse samples, validated for flow cytometry. Ideal for cell adhesion and receptor signaling studies using FC.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Mouse Integrin alpha 5 (ITGA5/CD49e) mAb
Remarks/Alias	FNRA, CD49e, VLA-5, VLA5A
Species	Human
GeneID (Human)	3678
GeneID	3678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-995 of human Integrin alpha 5 (ITGA5/CD49e) (NP_002196.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	P08648
Protein Weight	115kDa
Shipping	Ice bag