

Anti-Integrin Alpha 5 (Itga5/Cd49E) Polyclonal Antibody For Wb And Elisa - P08648

Item Number: CM0021116



Introduction

Rabbit polyclonal antibody against human Integrin alpha 5 (ITGA5/CD49e) validated for Western blot and ELISA. Cross-reacts with mouse. Suitable for studying fibronectin receptor function in cell adhesion and signaling.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 5 (ITGA5/CD49e) Rabbit pAb
Remarks/Alias	FNRA; CD49e; VLA-5; VLA5A; Integrin alpha 5 (ITGA5/CD49e)
Species	Human
GeneID (Human)	3678
GeneID	3678
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08648
Protein Weight	115kDa
Shipping	Ice bag