

Anti-Integrin Alpha 5 (Itga5/Cd49E) Rabbit Mab For Wb, Ihc-P, Fc (Intra), Elisa - P08648

Item Number: CM0002023



Introduction

Rabbit monoclonal ITGA5 antibody validated for WB, IHC-P, flow cytometry (intracellular) and ELISA. Cross-reacts with human, mouse and rat integrin alpha-5 fibronectin receptor subunit.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 5 (ITGA5/CD49e) Rabbit mAb
Remarks/Alias	FNRA; CD49e; VLA-5; VLA5A; Integrin alpha 5 (ITGA5/CD49e)
Species	Human
GeneID (Human)	3678
GeneID	3678
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 950-1049 of human Integrin alpha 5 (ITGA5/CD49e) (P08648).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08648
Protein Weight	115kDa
Shipping	Ice bag