

Anti-Integrin Alpha 6 (Itga6/Cd49F) Rabbit Polyclonal Antibody For Wb And Elisa - P23229

Item Number: CM0028388



Introduction

Rabbit polyclonal antibody against human Integrin alpha 6 (ITGA6/CD49f), validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Suitable for laminin receptor, cell adhesion, and hemidesmosome research.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 6 (ITGA6/CD49f) Rabbit pAb
Remarks/Alias	JEB6; CD49f; VLA-6; ITGA6A; ITGA6B; Integrin alpha 6 (ITGA6/CD49f)
Species	Human
GeneID (Human)	3655
GeneID	3655
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 481-591 of human Integrin alpha 6 (ITGA6/CD49f) (NP_001073286.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23229
Protein Weight	127kDa
Shipping	Ice bag