

Anti-Integrin Alpha-6/Itga6 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P23229

Item Number: CM0005011



Introduction

Recombinant rabbit monoclonal antibody against human Integrin alpha-6/ITGA6 (CD49f, VLA-6). Suitable for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for laminin receptor, cell adhesion, and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	Anti-Integrin alpha-6/ITGA6 Rabbit mAb
Remarks/Alias	JEB6; CD49f; VLA-6; ITGA6A; ITGA6B
Species	Human
Gene ID (Human)	3655
Gene ID	3655
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 505-663 of human Integrin alpha-6/ITGA6 (NP_000201.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23229
Protein Weight	127kDa
Shipping	Ice bag