

Anti-Integrin Alpha-V (Cd51) Monoclonal Antibody For Flow Cytometry (Intracellular) - P06756

Item Number: CM0008388



Introduction

-conjugated rabbit anti-human CD51 (integrin alpha-V) monoclonal antibody for intracellular flow cytometry. Validated in human samples, it detects a 116kDa receptor involved in cell adhesion, migration, and viral entry via RGD ligand binding.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD51 mAb
Remarks/Alias	CD51; MSK8; VNRA; VTNR
Species	Human
GeneID	3685
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 949-1048 of human CD51 (NP_002201.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P06756
Protein Weight	116kDa
Shipping	Ice bag