

Anti-Integrin Alpha V (Itgav/Cd51) Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P06756

Item Number: CM0022335



Introduction

Rabbit polyclonal antibody against human Integrin alpha V (ITGAV/CD51), validated for WB, IHC-P, IP, and ELISA; recognizes human, mouse, and rat. Ideal for studying cell adhesion, integrin signaling, and viral entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha V (ITGAV/CD51) Rabbit pAb
Remarks/Alias	CD51; MSK8; VNRA; VTNR; Integrin alpha V (ITGAV/CD51)
Species	Human
GeneID (Human)	3685
GeneID	3685
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06756
Protein Weight	116kDa
Shipping	Ice bag