

Anti-Integrin Beta 1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09055

Item Number: CM0027168



Introduction

Rabbit polyclonal antibody against mouse Integrin beta 1 (CD29), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human. Ideal for cell adhesion and integrin signaling research.

[Learn More](#)

Parameter	Value
Product Name	Integrin beta 1 Rabbit pAb
Remarks/Alias	CD29, Fnrb, gpIIa, Gm9863, 4633401G24Rik
Species	Mouse
Gene ID (Human)	3688
Gene ID	16412
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-728 of mouse Integrin beta 1 (NP_034708.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09055
Protein Weight	88kDa
Shipping	Ice bag