

Anti-Integrin Beta 3 (Itgb3/Cd61) Rabbit Monoclonal Antibody For Wb, Ip, Fc, Elisa - P05106

Item Number: CM0003706



Introduction

High-quality rabbit monoclonal antibody against human Integrin beta 3 (ITGB3/CD61), validated for WB, IP, FC (intra), and ELISA. Recognizes human, mouse, and rat ITGB3; useful for platelet adhesion and integrin signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Integrin beta 3 (ITGB3/CD61) Rabbit mAb
Remarks/Alias	GT, GT2, CD61, GP3A, BDPLT2, GPIIIa, BDPLT16, BDPLT24, Integrin beta 3 (ITGB3/CD61)
Species	Human
GeneID (Human)	3690
GeneID	3690
Immunogen	A synthetic peptide corresponding to amino acids 689-788 of human Integrin beta 3 (P05106)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, FC (intra), ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05106
Protein Weight	87kDa
Shipping	Ice bag