

Anti-Integrin Beta 3 (Itgb3/Cd61) Rabbit Polyclonal Antibody For Wb, Elisa - P05106

Item Number: CM0018343



Introduction

Integrin beta 3 (ITGB3/CD61) rabbit polyclonal antibody for WB and ELISA. Recognizes human and mouse ITGB3, mediating platelet aggregation and integrin-mediated adhesion. Recombinant protein immunogen.

[Learn More](#)

Parameter	Value
Product Name	Integrin beta 3 (ITGB3/CD61) Rabbit pAb
Remarks/Alias	GT; GT2; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16; BDPLT24; Integrin beta 3 (ITGB3/CD61)
Species	Human
GeneID (Human)	3690
GeneID	3690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 609-788 of human ITGB3 (NP_000203.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05106
Protein Weight	87kDa
Shipping	Ice bag