

Anti-Integrin Beta 3 (Itgb3/Cd61) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05106

Item Number: CM0021477



Introduction

Rabbit polyclonal antibody against integrin beta 3 (ITGB3/CD61), validated for WB, IHC-P, and ELISA with human, mouse, rat reactivity. Perfect for platelet aggregation and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Integrin beta 3 (ITGB3/CD61) Rabbit pAb
Remarks/Alias	GT; GT2; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16; BDPLT24; Integrin beta 3 (ITGB3/CD61)
Species	Human
GeneID (human)	3690
GeneID	3690
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 80-430 of human Integrin beta 3 (ITGB3/CD61) (NP_000203.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05106
Protein Weight	87kDa
Shipping	Ice bag