

Anti-Integrin Beta 4 (Itgb4/Cd104) Monoclonal Antibody For If-P, Ihc-P, Elisa - P16144

Item Number: CM0007436



Introduction

Rabbit monoclonal antibody against human Integrin beta-4 (ITGB4/CD104), validated for IF-P, IHC-P, ELISA, with cross-reactivity to mouse and rat. This laminin receptor is critical for epithelial cell adhesion and hemidesmosome integrity.

[Learn More](#)

Parameter	Value
Product Name	Integrin beta 4 (ITGB4/CD104) Rabbit mAb
Remarks/Alias	CD104; GP150; JEB5A; JEB5B; Integrin β 4
Species	Human
Gene ID (Human)	3691
Gene ID	3691
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16144
Protein Weight	202kDa
Shipping	Ice bag