

Anti-Integrin Beta 4 (Itgb4/Cd104) Monoclonal Antibody For Wb, Elisa - P16144

Item Number: CM0011051



Introduction

Recombinant rabbit monoclonal antibody targeting human Integrin beta 4 (ITGB4/CD104). Suitable for WB and ELISA; cross-reacts with mouse. Critical for hemidesmosome integrity and keratinocyte polarity.

[Learn More](#)

Parameter	Value
Product Name	Integrin beta 4 (ITGB4/CD104) Rabbit mAb
Remarks/Alias	CD104; GP150; JEB5A; JEB5B; Integrin beta 4 (ITGB4/CD104)
Species	Human
Gene ID	3691
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1604-1822 of human ITGB4 (P16144)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16144
Protein Weight	202kDa
Shipping	Ice bag