


Anti-Integrin Beta 4 (Itgb4/Cd104) Rabbit Polyclonal Antibody For Wb, Elisa - P16144

Item Number: CM0015433



Introduction

Rabbit polyclonal antibody against human Integrin beta 4 (ITGB4/CD104), validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Suitable for hemidesmosome and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	Integrin beta 4 (ITGB4/CD104) Rabbit pAb
Remarks/Alias	CD104; GP150; JEB5A; JEB5B; Integrin beta 4 (ITGB4/CD104)
Species	Human
Gene ID	3691
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 611-710 of human Integrin beta 4 (ITGB4/CD104) (NP_000204.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16144
Protein Weight	202kDa
Shipping	Ice bag