

Anti-Integrin Beta 8 Monoclonal Antibody For Wb, Elisa - P26012

Item Number: CM0005331

Introduction

High-quality rabbit monoclonal antibody targeting human Integrin beta 8 (ITGB8). Validated for Western blot and ELISA applications with cross-reactivity to human, mouse, and rat. Ideal for investigating integrin-mediated cell adhesion, signaling, and TGF-beta activation pathways.

[Learn More](#)

Parameter	Value
Product Name	Integrin beta 8 Rabbit mAb
Remarks/Alias	ITGB8
Species	Human
GeneID (Human)	3696
GeneID	3696
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-320 of human Integrin beta 8 (NP_002205.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26012
Protein Weight	71 kDa/86 kDa
Shipping	Ice bag