

Anti-Integrin-B1/Cd29 Monoclonal Antibody For If/Icc, Elisa - P05556

Item Number: CM0002122



Introduction

Recombinant rabbit monoclonal antibody targeting human Integrin- β 1/CD29 (P05556). Validated for IF/ICC and ELISA, with cross-reactivity to mouse and rat. Ideal for cell adhesion and integrin signaling research.

[Learn More](#)

Parameter	Value
Product Name	Integrin- β 1/CD29 Rabbit mAb
Remarks/Alias	CD29; FNRB; MDF2; VLAB; GPIIA; MSK12; VLA-BETA; Integrin- β 1/CD29
Species	Human
Gene ID (Human)	3688
Gene ID	3688
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-728 of human Integrin- β 1/CD29 (NP_002202.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05556
Protein Weight	88kDa
Shipping	Ice bag