

Anti-Integrin-B1/Cd29 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05556

Item Number: CM0020043



Introduction

High-quality rabbit polyclonal antibody against human Integrin-β1/CD29 (P05556). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for cell adhesion and integrin signaling research.

[Learn More](#)

Parameter	Value
Product Name	Integrin-β1/CD29 Rabbit pAb
Remarks/Alias	CD29; FNRB; MDF2; VLAB; GPIIA; MSK12; VLA-BETA; Integrin-β1/CD29
Species	Human
Gene ID (Human)	3688
Gene ID	3688
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-728 of human Integrin-β1/CD29 (NP_002202.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05556
Protein Weight	88kDa
Shipping	Ice bag