

Anti-Integrin-B1/Cd29 Rabbit Polymab® - P05556

Item Number: CM0003214



Introduction

Rabbit monoclonal antibody against human Integrin-β1/CD29 (P05556) validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human and mouse. Ideal for cell adhesion and viral entry studies.

[Learn More](#)

Parameter	Value
Product Name	Integrin-β1/CD29 Rabbit PolymAb®
Remarks/Alias	CD29; FNRB; MDF2; VLAB; GPIIA; MSK12; VLA-BETA
Species	Human
Gene ID (Human)	3688
Gene ID	3688
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 1-728 of human Integrin-β1/CD29 (NP_002202.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05556
Protein Weight	88kDa
Shipping	Ice bag