


Anti-Integrin B5/Itgb5 Monoclonal Antibody For Fc - P18084

Item Number: CM0005555



Introduction

Anti-ITGB5 rabbit monoclonal antibody specifically recognized human Integrin beta-5. Validated for flow cytometry (FC), it targets the extracellular domain for cell surface detection. Research use in adhesion and viral entry studies.

[Learn More](#)

Parameter	Value
Product Name	Integrin β 5/ITGB5 Rabbit mAb
Remarks/Alias	ITGB5; integrin beta-5; Integrin β 5/ITGB5
Species	Human
GeneID (Human)	3693
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-490 of human Integrin β 5/ITGB5 (NP_002204.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P18084
Protein Weight	88kDa
Shipping	Ice bag