


Anti-Ints10 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nvr2

Item Number: CM0026917



Introduction

High-quality rabbit polyclonal antibody against human Integrator complex subunit 10 (INTS10), validated for western blot and ELISA. Reacts with human, mouse, and rat samples, ideal for transcription termination research.

[Learn More](#)

Parameter	Value
Product Name	INTS10 Rabbit pAb
Remarks/Alias	INT10; C8orf35; INTS10
Species	Human
GeneID (Human)	55174
GeneID	55174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 451-710 of human INTS10 (NP_060612.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVR2
Protein Weight	82kDa
Shipping	Ice bag