


Anti-Ints11 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Ta45

Item Number: CM0023190



Introduction

High-quality rabbit polyclonal antibody targeting human INTS11 (Integrator complex subunit 11), suitable for Western blot and ELISA applications. Shows cross-reactivity with human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	INTS11 Rabbit pAb
Remarks/Alias	RC68; INT11; RC-68; CPSF3L; CPSF73L; INTS11
Species	Human
GeneID (human)	54973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 331-600 of human INTS11 (NP_060341.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5TA45
Protein Weight	68kDa
Shipping	Ice bag