

Anti-Ints4 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q96Hw7

Item Number: CM0031061



Introduction

High-quality rabbit polyclonal antibody targeting human INTS4 (Integrator complex subunit 4) for WB, IP, ELISA applications. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	INTS4 Rabbit pAb
Remarks/Alias	INT4, MST093, INTS4
Species	Human
GeneID	92105
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human INTS4 (NP_291025.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96HW7
Protein Weight	108kDa
Shipping	Ice bag