


Anti-Ints5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6P9B9

Item Number: CM0013227



Introduction

Rabbit polyclonal antibody against human INTS5 (Q6P9B9) validated for WB, IHC-P, ELISA. Reacts with human and mouse integrator complex subunit 5, a key regulator of RNA polymerase II transcription termination.

[Learn More](#)

Parameter	Value
Product Name	INTS5 Rabbit pAb
Remarks/Alias	INT5; KIAA1698; INTS5
Species	Human
Gene ID (Human)	80789
Gene ID	80789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-400 of human INTS5 (NP_085131.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6P9B9
Protein Weight	108kDa
Shipping	Ice bag