


## Anti-Ints6 Polyclonal Antibody For Wb, Elisa - Q9UI03

Item Number: CM0029454



### Introduction

Rabbit polyclonal antibody against human INTS6 (Integrator complex subunit 6), validated for WB and ELISA. Cross-reacts with mouse INTS6. Useful for transcription termination and tumor suppression studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | INTS6 Rabbit pAb                                                                                                                     |
| Remarks/Alias    | HDB; INT6; DBI-1; DDX26; DICE1; DDX26A; Notch12; INTS6                                                                               |
| Species          | Human                                                                                                                                |
| Gene ID (Human)  | 26512                                                                                                                                |
| Gene ID          | 26512                                                                                                                                |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human INTS6 (NP_036273.1) |
| Source           | Rabbit                                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                                |
| Application      | WB, ELISA                                                                                                                            |
| Cross Reactivity | Human, Mouse                                                                                                                         |
| SWISS            | Q9UI03                                                                                                                               |
| Protein Weight   | 100kDa                                                                                                                               |
| Shipping         | Ice bag                                                                                                                              |