


Anti-Ints6 Rabbit Polyclonal Antibody For Wb, Elisa - Q9UI03

Item Number: CM0022919



Introduction

Rabbit polyclonal antibody targeting human INTS6, an integrator complex subunit involved in transcription termination. Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for cancer and transcription research.

[Learn More](#)

Parameter	Value
Product Name	INTS6 Rabbit pAb
Remarks/Alias	HDB; INT6; DBI-1; DDX26; DICE1; DDX26A; Notch12; INTS6
Species	Human
GeneID (Human)	26512
GeneID	26512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-887 of human INTS6 (NP_036273.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UI03
Protein Weight	100kDa
Shipping	Ice bag