


Anti-Ints7 Polyclonal Antibody For Wb, Elisa - Q9Nvh2

Item Number: CM0024616



Introduction

Rabbit polyclonal antibody against human INTS7, validated for Western blot and ELISA. Reacts with human and mouse integrator complex subunit 7, a transcriptional termination regulator.

[Learn More](#)

Parameter	Value
Product Name	INTS7 Rabbit pAb
Remarks/Alias	INT7; C1orf73; INTS7
Species	Human
GeneID (Human)	25896
GeneID	25896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 702-948 of human INTS7 (NP_001186740.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NVH2
Protein Weight	107kDa
Shipping	Ice bag