


Anti-Ints8 Rabbit Polyclonal Antibody For Wb, Elisa - Q75Qn2

Item Number: CM0022484



Introduction

Anti-INTS8 rabbit polyclonal antibody, validated for WB and ELISA, reacts with human and mouse samples. Ideal for research on integrator complex-mediated transcription termination and RNA processing.

[Learn More](#)

Parameter	Value
Product Name	INTS8 Rabbit pAb
Remarks/Alias	INT8; NEDCHS; C8orf52; INTS8
Species	Human
GeneID (Human)	55656
GeneID	55656
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-800 of human INTS8 (NP_060334.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q75QN2
Protein Weight	113kDa
Shipping	Ice bag