

Anti-Involucrin Rabbit Monoclonal Antibody For Wb - P07476

Item Number: CM0003710

Introduction

Recombinant rabbit monoclonal antibody targeting human Involucrin. Validated for WB, IHC-P, and ELISA; cross-reacts with rat. Suited for keratinocyte differentiation and cornified envelope studies.

[Learn More](#)

Parameter	Specification
Product Name	Involucrin Rabbit mAb
Species	Human
Gene ID	3713
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human Involucrin (P07476)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P07476
Protein Weight	68 kDa
Shipping	Ice bag