

Anti-Invs Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y283

Item Number: CM0016307

Introduction

Rabbit polyclonal antibody against human INVS (Inversin), validated for Western blotting and ELISA. Cross-reacts with human and mouse. Suitable for studies on renal development, left-right axis determination, and Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	INVS Rabbit pAb
Remarks/Alias	INV; NPH2; NPHP2; INVS
Species	Human
GeneID (Human)	27130
GeneID	27130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human INVS (NP_055240.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y283
Protein Weight	118kDa
Shipping	Ice bag