

Anti-IP3R Rabbit Polyclonal Antibody For Wb, Elisa - Q14643

Item Number: CM0025529

Introduction

High-quality rabbit polyclonal antibody specific for IP3R (ITPR1). Validated for WB and ELISA applications, with cross-reactivity to human and mouse. Ideal for investigating calcium signaling, apoptosis, and channel function.

[Learn More](#)

Parameter	Value
Product Name	IP3R Rabbit pAb
Remarks/Alias	ACV; CLA4; IP3R; IP3R1; SCA15; SCA16; SCA29; INSP3R1; PPP1R94
Species	Human
GeneID (human)	3708
GeneID	3708
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1812-1911 of human IP3R (NP_001365381).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14643
Protein Weight	314kDa
Shipping	Ice bag