

Anti-IP6K1 Polyclonal Antibody For Wb, IHC-P, Elisa - Q92551

Item Number: CM0014347

Introduction

Rabbit polyclonal antibody raised against recombinant human IP6K1 protein. Validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying inositol phosphate signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	IP6K1 Rabbit pAb
Remarks/Alias	PIUS; IHPK1; IP6K1
Species	Human
Gene ID (Human)	9807
Gene ID	9807
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human IP6K1 (NP_695005.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92551
Protein Weight	50kDa
Shipping	Ice bag