

Anti-Ip6K1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92551

Item Number: CM0004455

Introduction

High-quality anti-IP6K1 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat IP6K1 (Q92551), ~50 kDa. Ideal for inositol phosphate metabolism research.

[Learn More](#)

Parameter	Value
Product Name	IP6K1 Rabbit mAb
Remarks/Alias	PiUS; IHPK1; IP6K1
Species	Human
Gene ID (Human)	9807
Gene ID	9807
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human IP6K1 (NP_695005.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92551
Protein Weight	50 kDa
Shipping	Ice bag