

Anti-IP6K2 Polyclonal Antibody For If/Icc, Elisa - Q9Uhh9

Item Number: CM0016974

Introduction

Rabbit polyclonal antibody targeting human IP6K2, suited for IF/ICC and ELISA. Useful for exploring inositol phosphate metabolism, InsP6 to InsP7 conversion, and nuclear signaling pathways.

[Learn More](#)

Parameter	Specification
Product Name	IP6K2 Rabbit pAb
Remarks/Alias	PIUS; IHPK2; InsP6K2; IP6K2
Species	Human
GeneID (Human)	51447
GeneID	51447
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human IP6K2 (NP_057375.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9UHH9
Protein Weight	49kDa
Shipping	Ice bag