

Anti-Ipa1 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Exz2

Item Number: CM0021231

Introduction

Rabbit polyclonal antibody against rice IPA1/SPL14, validated for WB and ELISA. Detects Oryza sativa IPA1 (~42 kDa). Ideal for plant architecture and defense research.

[Learn More](#)

Parameter	Value
Product Name	IPA1 Rabbit pAb
Remarks/Alias	WFP, IPA1, SPL14, OsSPL14
Species	Oryza sativa
GeneID	4345998
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 209-417 of Oryza sativa IPA1 (Q7EXZ2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q7EXZ2
Protein Weight	42kDa
Shipping	Ice bag