


Anti-Ipo5 Rabbit Polyclonal Antibody For Wb, Elisa - 000410

Item Number: CM0021592



Introduction

Anti-IPO5 rabbit polyclonal antibody for WB and ELISA applications. Detects importin-5 in human, mouse, and rat samples. Ideal for studying nuclear transport and nuclear localization signals.

[Learn More](#)

Parameter	Value
Product Name	IPO5 Rabbit pAb
Remarks/Alias	IMB3; Pse1; imp5; KPNB3; RANBP5; IPO5
Species	Human
GeneID (human)	3843
GeneID	3843
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human IPO5 (NP_002262.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00410
Protein Weight	124kDa
Shipping	Ice bag