

Anti-Ipo7 Polyclonal Antibody For Wb, If/Icc, Elisa - 095373

Item Number: CM0021819

Introduction

Rabbit polyclonal antibody targeting human importin-7 (IPO7/RanBP7). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat samples. Ideal for nuclear transport and importin pathway research.

[Learn More](#)

Parameter	Value
Product Name	IPO7 Rabbit pAb
Remarks/Alias	Imp7; RANBP7; IPO7
Species	Human
Gene ID (Human)	10527
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 900 to the C-terminus of human IPO7 (NP_006382.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95373
Protein Weight	120kDa
Shipping	Ice bag