

Anti-Ipo8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15397

Item Number: CM0018450

Introduction

High-quality IPO8 Rabbit pAb (alias RANBP8) validated for Western blot, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for nuclear transport and importin research.

[Learn More](#)

Parameter	Value
Product Name	IPO8 Rabbit pAb
Remarks/Alias	VISS; RANBP8; IPO8
Species	Human
Gene ID (Human)	10526
Gene ID	10526
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 878-1037 of human IPO8 (NP_006381.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15397
Protein Weight	120kDa
Shipping	Ice bag