


Anti-Ipo9 Rabbit Polyclonal Antibody For Wb And Elisa - Q96P70

Item Number: CM0013926



Introduction

High-quality rabbit polyclonal antibody against human IPO9 (Importin-9). Validated for WB and ELISA applications, with cross-reactivity to human and mouse. Ideal for nucleocytoplasmic transport and histone chaperone studies.

[Learn More](#)

Parameter	Value
Product Name	IPO9 Rabbit pAb
Remarks/Alias	Imp9; IPO9
Species	Human
GeneID (Human)	55705
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 827-1041 of human IPO9 (NP_060555.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96P70
Protein Weight	116kDa
Shipping	Ice bag