

Anti-Iqgap1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P46940

Item Number: CM0010167



Introduction

IQGAP1 Rabbit mAb for IVD assay development: validated for WB, IHC-P, ELISA with human/mouse reactivity. Targets Ras GTPase-activating-like protein IQGAP1 (p195, SAR1), involved in actin dynamics and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	IQGAP1 Rabbit mAb
Remarks/Alias	SAR1; p195; HUMORFA01; IQGAP1
Species	Human
GeneID (human)	8826
GeneID	8826
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human IQGAP1 (NP_003861.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46940
Protein Weight	189kDa
Shipping	Ice bag