

Anti-Iqgap1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P46940

Item Number: CM0028903



Introduction

CamelBio offers anti-IQGAP1 rabbit polyclonal antibody (P46940) validated for WB, IHC-P, IF/ICC, IP, ELISA in human, mouse, rat. Ideal for studying actin cytoskeleton dynamics and cell signaling.

[Learn More](#)

Parameter	Value
Product Name	IQGAP1 Rabbit pAb
Remarks/Alias	SAR1; p195; HUMORFA01; IQGAP1
Species	Human
GeneID (Human)	8826
GeneID	8826
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 500-750 of human IQGAP1 (NP_003861.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46940
Protein Weight	189kDa
Shipping	Ice bag