

Anti-Iqgap2 Monoclonal Antibody For Wb - Q13576

Item Number: CM0003864



Introduction

High-quality Anti-IQGAP2 Rabbit Monoclonal Antibody for WB, IHC-P, ELISA. Recognizes human, mouse, rat IQGAP2, a 181 kDa Ras GTPase-activating-like protein binding CDC42/RAC1 and calmodulin.

[Learn More](#)

Parameter	Value
Product Name	IQGAP2 Rabbit mAb
Remarks/Alias	IQGAP2
Species	Human
Gene ID (Human)	10788
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1301-1495 of human IQGAP2 (Q13576).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13576
Protein Weight	181kDa
Shipping	Ice bag