

Anti-Iqgap2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13576

Item Number: CM0016150



Introduction

Rabbit polyclonal antibody targeting human IQGAP2, validated for WB, IHC-P, and ELISA applications. Reacts with human, mouse, and rat samples. Useful for studying GTPase signaling and calmodulin interactions.

[Learn More](#)

Parameter	Value
Product Name	IQGAP2 Rabbit pAb
Remarks/Alias	IQGAP2
Species	Human
GeneID (Human)	10788
GeneID	10788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1330-1440 of human IQGAP2 (NP_006624.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13576
Protein Weight	181kDa
Shipping	Ice bag