


Anti-Irak1 Polyclonal Antibody For Wb, If/Icc, Elisa - P51617

Item Number: CM0017520



Introduction

High-quality rabbit polyclonal antibody targeting human IRAK1, essential for innate immune response research via TLR/IL-1R pathways. Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for diagnostics and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	IRAK1 Rabbit pAb
Remarks/Alias	IRAK; pelle; IRAK1
Species	Human
GeneID (human)	3654
GeneID	3654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 424-633 of human IRAK1 (NP_001020414.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51617
Protein Weight	77kDa
Shipping	Ice bag