


Anti-Irak1 Rabbit Monoclonal Antibody For Wb, Elisa - P51617

Item Number: CM0008173



Introduction

High-quality anti-IRAK1 rabbit monoclonal antibody validated for Western blot and ELISA. Targets human IRAK1 (P51617), a key innate immunity kinase involved in TLR/IL-1R signaling. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	IRAK1 Rabbit mAb
Remarks/Alias	IRAK; pelle; IRAK1
Species	Human
GeneID (Human)	3654
GeneID	3654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 215-499 of human IRAK1 (P51617)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P51617
Protein Weight	77kDa
Shipping	Ice bag