

Anti-Irak2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O43187

Item Number: CM0004110

Introduction

High-performance rabbit monoclonal antibody targeting human IRAK2, validated for WB, IHC-P, and ELISA. Detects human and rat IRAK2, a mediator of IL-1 receptor signaling. Ideal for inflammation and immunology research.

[Learn More](#)

Parameter	Value
Product Name	IRAK2 Rabbit mAb
Remarks/Alias	IRAK-2; IRAK2
Species	Human
Gene ID (Human)	3656
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-250 of human IRAK2 (O43187)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O43187
Protein Weight	69kDa
Shipping	Ice bag