


Anti-Irak3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y616

Item Number: CM0028819



Introduction

IRAK3 Rabbit Polyclonal Antibody detects IRAK3 (IRAK-M) in human, mouse, rat samples. Suitable for WB, IF/ICC, ELISA. IRAK3 is an inactive kinase that regulates TLR/IL-1R signaling, relevant in inflammation and asthma research.

[Learn More](#)

Parameter	Value
Product Name	IRAK3 Rabbit pAb
Remarks/Alias	ASRT5, IRAKM, IRAK3
Species	Human
GeneID	11213
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 367-596 of human IRAK3 (NP_009130.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y616
Protein Weight	68kDa
Shipping	Ice bag