


Anti-Irak4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nwz3

Item Number: CM0024586



Introduction

Rabbit polyclonal antibody against human IRAK4, a key kinase in innate immunity. Validated for WB, IF/ICC, and ELISA; reacts with human, mouse, and rat samples. Ideal for TLR/IL-1R signaling research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	IRAK4 Rabbit pAb
Remarks/Alias	IPD1; IMD67; REN64; IRAK-4; NY-REN-64; IRAK4
Species	Human
GeneID (human)	51135
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWZ3
Protein Weight	52kDa
Shipping	Ice bag