


Anti-Irak4 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nwz3

Item Number: CM0002419



Introduction

KD-validated recombinant rabbit monoclonal antibody against human IRAK4 for Western blot and ELISA. Suitable for studying IRAK4's role in Myddosome complex formation, TLR/IL-1R signaling, and NF-kB activation.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] IRAK4 Rabbit mAb
Remarks/Alias	IPD1; IMD67; REN64; IRAK-4; NY-REN-64; [KD Validated] IRAK4
Species	Human
GeneID (Human)	51135
GeneID	51135
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NWZ3
Protein Weight	52kDa
Shipping	Ice bag