


Anti-Irak4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nwz3

Item Number: CM0026430



Introduction

IRAK4 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA, recognizing human, mouse, and rat IRAK4. Ideal for innate immune response and TLR/IL-1R signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	IRAK4 Rabbit pAb
Remarks/Alias	IPD1; IMD67; REN64; IRAK-4; NY-REN-64; IRAK4
Species	Human
Gene ID	51135 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human IRAK4 (NP_057207.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWZ3
Protein Weight	52kDa
Shipping	Ice bag